



All India Institute of Medical Sciences, Bhubaneswar
At - Sijua (Patrapada), Post - Dumuduma, Bhubaneswar (Odisha) - 751 019

Reference No : J-11028(051)/2017/S&P /55

Dated : 14th Sep 2017

**Notice inviting Supply of Semi Motorized Rotary Microtome
as proprietary Items at AIIMS, Bhubaneswar**

Name of the Equipment : **Semi Motorized Rotary Microtome, (Mfg. by Lieca
Mikrosysteme Vertrieb GmbH, Germany).**

All India Institute of Medical Sciences (AIIMS), Bhubaneswar, Odisha, invites comments regarding supply of Consumables/ Accessories for **Semi Motorized Rotary Microtome** under Proprietary Article Certified by **M/s Lieca Mikrosysteme Vertrieb GmbH**, Germany. The proposal submitted by the firm (M/s DHR Holding India Pvt. Ltd., Kolkata) and Proprietary Article Certified, Authorized Certificate and quoted price and Specification are attached & uploaded on website for comments.

Reputed firms, Supplier and Manufacturer are requested to submit their Comments or best offer along with complete details of specifications for **Semi Motorized Rotary Microtome**. The comments or any information should reach the Institute on or before **23-09-2017 at 12.00 PM** falling which it will be presumed that no other vendor is having any comment to offer and case will be decided on merit.

**Sr. Procurement-cum Store Officer
AIIMS, BBSR**

Copy to:

- | | | |
|----------------------------|---|------------------------------|
| 1. HOD, Dept. of Pathology | : | for information Please. |
| 2. Accounts Officer | : | for information Please. |
| 3. IT Cell | : | for hoisting in the website. |

Living up to Life



Leica Mikrosysteme Vertrieb GmbH
Ernst Leitz Strasse 17-37
35578 Wetzlar
Germany

Date: 11.09.2017

To,

The Director

All India Institute of Medical Sciences, Bhubaneswar
Sijua, Near Bijupatnayak Police Academy,
Patrapada, Bhubaneswar, Odisha 751019

Dear Sir/Madam

It is hereby declared and confirmed that Leica Semi Motorized Rotary Microtome Model RM 2245, vide Quotation No. SB0014 is our proprietary item and to the best of our knowledge no company other than us manufacture this product anywhere else in the world.

For LEICA MIKROSYSTEME VERTRIEB GmbH.

Authorized Signatory

Qtn No.:SB 0014/12092017

Date: 12.09.2017

To,

Sub : **Special Discounted Offer for Leica Instrument**

Leica Manual Microtome-Model RM 2245

Sl No.	Description
1.	• Leica RM2245 Semi Motorized Microtome Designed for routine and research applications in histopathology and industrial quality assurance laboratories Leica RM2245 - Configuration 2 High-profile disposable blades 818

FOR destination Price	INR	1650000.00
Special Discount	INR	400000.00
Special Discounted FOR price	INR	1250000.00

Living up to Life

Leica
MICROSYSTEMS

19

Terms & Conditions (for Purchase in INR)

Currency: Indian Rupees (INR)

Terms: F.O.R. destination Price

Tax: GST @ 18% applicable extra

Payment Terms: As per AIIMS Bhubaneswar standard payment terms.

Delivery: Within 6-8 weeks from the date of receipt of your P.O.


Warranty: All the above quoted instrument will carry comprehensive warranty for 5 years from the date of installation.

Installation: Free of cost at your site.

Offer Validity: 60 days.

Order on: Please note order and supply will be made by our authorized dealer M/S Bioapps.

For DHR Holding India Pvt. Ltd.



Sinchan Banerjee
Area Sales Manager
Leica Biosystems



SI No	Unique Feature	Description	Customer Specific Benefit
1	Fixture for specimen clamps - precision orienting	8 degree defined X/Y fine adjustment orientation with click stops every 2 degree and clear zero position with red indicators	Help in processing High sample load with precision
2	Patented adjustable handwheel balance	1. Smoothly running handwheel when using manually. 2. No need of 2 kg heavy balance weight in handwheel.	1. Reduced risk of ergonomic injury when working manually with high sample load 2. Increased specimen safety as object does not drop into the knife when stopping sectioning
		3. Ensures always optimal stopping point of object head.	
		3. Possibility to balance heavy clamps - such as Supermegacassette	Usefull for processing big tissue samples.
3	Selection between: continuous and step feeding in sectioning mode	1. Continuous feeding: same function as in trimming mode. Fast movement back and forth.	Continuous feeding: no confusion of the reaction of the coarse feed butting in trimming and sectioning mode. The coarse feed could however be too crude for the sectioning mode.
		2. Step feeding: After trimming only little steps are needed to bring the specimen back to the knife.	2. Step feeding: For improved specimen safety: This mode is recommended for very small and sensitive specimens which could easily be damaged by too crude course feeding.

4	Rocking mode	No need to co-ordinate handwheel up and down movement and simultaneously pressing coarse feed buttons.	Avoids stress. Reduce strain on the technician and increase productivity.
5	Large ergonomic handwheel handle	ergonomically shaped handwheel handle for firm grip of the whole hand	Firm grip - avoids ergonomic stress
6	Ruler on microtome base	1. Ruler helps to find ideal knife holder position for different mold sizes.	Time saving and easy to use
		2. Use together with object exchange in home position	

Feature	Description	Benefit
<u>Knife holder NZ</u>	Stable central clamping pressure plate allows usage of the entire length of expensive c or d profile reusable tungsten carbide knives.	High section quality of hard specimens
		Economics: Entire length of the knife can be used
Vertical Stroke 70 mm :	Large vertical stroke	1. Allows sectioning of Supermegacassettes. 2. Ensures safe changing of specimens as a result of a great distance to the knife edge compared to competitor microtomes
Trim thickness range : 1 to 600um	wide range	Users can select preferred value for their application
Section thickness range : 0,5 - 100um	wide range	

	Semi-thin sectioning up to 0.5um	Users can select preferred value for their application
<u>Fixed knife or blade holder base</u>	High stability - no integrated lateral displacement function weakening the system	High section quality,
	Cut-out for back-lighting system	
Coarse feed buttons	4 conveniently positioned coarse feed buttons on separate control panel.	Intuitive operation
	2 speeds per direction.	
	Buttons can be easily found due to touch sensitive and color marked buttons.	

From Eye to Insight



Leica Mikrosysteme Vertrieb GmbH · Ernst-Leitz-Str. 17 · 37 · 35578 Wetzlar, Germany

To whom it may concern

Wetzlar, February 2017

Distribution Authorization

Since November 2011 Leica Microsystems is operating directly in India through a legal entity of its parent company Danaher Corporation called M/S DHR Holding India Pvt. Ltd.

Therefore, by this document we confirm that:

DHR Holding India Pvt. Ltd.
3rd Floor, B Wing,
Art Guild House,
Phoenix Market City, LBS Road,
Kurla West, Mumbai 400070
INDIA

has the exclusive authorization to distribute our products of Compound Microscopy, Stereo Microscopy, Confocal Microscopy, Surgical Microscopy and EM Sample Preparation product lines manufactured by Leica Microsystems.

and is our exclusively authorized distributor for the Histology Systems, Histology Consumables, Advanced Staining Systems, Pathology Imaging: Clinical Microscopy, Digital pathology & Cytovision product lines manufactured by Leica Biosystems.

As such, DHR Holding India Pvt. Ltd. is authorized to offer for sales and servicing all said products throughout the territory of India, Nepal, Bhutan, Bangladesh & Sri Lanka.

This authorisation letter is valid until February 28, 2018 and may be prolonged in writing.

Leica Mikrosysteme Vertrieb GmbH, DSA

Markus Vornhagen
Leica Group Controller

Karim Sarace
Regional Manager Shared Services

Leica Mikrosysteme Vertrieb GmbH · GmbH mit Sitz in Wetzlar · Amtsgericht Wetzlar HRB 5464 ·

Geschäftsführer: Frank Bunge · Melanie Hammerschmidt-Broman · Markus Vornhagen

www.leica-microsystems.com

Deutsche Bank AG, Frankfurt · BIC DEUTDEFFXXX · IBAN DE34 5007 0010 0040 5100 00

Living up to Life

Leica

MICROSYSTEMS

DHR Holding India Pvt Ltd,
Arcadia, 3rd & 4th Floor, Premise No.31 Dr.Ambedkar Sarani, Kolkata, West Bengal - 700046
Telephone: 033 6637 6148/6150
Fax: 033 6637 6665
Email: contact.india@leica-microsystems.com

Ref: CL-17/Bioap/BBSR/005

Date: 11-09-2017

WHOM SO EVER IT MAY CONCERN

Subject: Authorization Letter

Dear Sir,

We, M/s.DHR Holding India Pvt. Ltd., being the exclusive Indian distributor for M/s. Leica Microsystems Vertrieb GmbH, Germany, hereby authorize:-

**M/S BIOAPPS
181, Bangur Avenue, Block-C
Kolkata-700055**

to offer, quote, negotiate and finalize the deal behalf of **Leica Biosystems, DHR Holding India Pvt. Ltd.**
Arcadia, 3rd & 4th Floor, Premise No.31, Dr. Ambedkar Sarani, Kolkata 700046, India

We assure you to take care of installation and after-sales-service for the products quoted and supplies by
M/S BIOAPPS, 181, Bangur Avenue, Block-C, Kolkata-700055

We are enclosing herewith the authorization certificate in our favor received from Leica for your ready reference.

Yours truly,
For, **DHR Holding India Pvt. Ltd.**


**Subhojit Dutta
Regional Sales Manager
Leica Biosystems Division.
Encl: As above.**



REGIONAL OFFICES AT:

Thane - 022 25986000, Delhi - 011 43766300, Lucknow - 0522 2346535, Chennai - 044 2347008, Hyderabad - 040 23313300,
Bangalore - 080 23230912, Thiruvananthapuram - 0471 2324064, Baroda - 0265 2342929.
www.leica-microsystems.com

Proprietary Article Certificate

Valid for the Current Financial Year

10

File No. and Date Reference :			
1	Description of article	Semi Motorized Rotary Microtome Model : RM 2245, Make : Leica Mikrosysteme Vertrieb GmbH	
2	Forecast of quantity /annual requirement	One time purchase	
3	Approximate estimated value for above quantity	Rs. 10,25,000.00	
4	Maker's (OEM) name and address	Leica Mikrosysteme Vertrieb GmbH Germany	
5	Name's of authorised dealers/ stockists		
6	I approve the above purchase on PAC basis and certify that : - Note- Tick to retain only one out of (b), C-1) or (c-2) whichever is applicable and cross out others. Please do confirm (a) by ticking it- without PAC certificate will be invalid.		
6 (a)	This is the only firm who is manufacturing / stocking this item. AND	Yes	
6 (b)	A Similar article is not manufactured / sold by any other firm, which could be used in lieu	No	
6 (c-1)	No other make/brand will be suitable for following tangible reasons (like OEM/Warranty spares.) : OR	No	
6 (c-2)	No other make/brand will be suitable and acceptable for following intangible reasons (if PAC was also given in the last procurement cycle, please also bring out efforts made since then to locate more sources) : <u>Justification is attached in Annexure-"A"</u>		
7	Reference of concurrence of finance wing to the proposal :		
History of PAC Purchase of this item for past three years may be given below : Not applicable			
Name of the Supplier			
Order/ Tender Reference & Date	Quantity Ordered	Basic Rate on Order (Rs.)	Adverse Performance Reported if Any

Smiti Randa Mishra
Signature of the Indenter/ User

Dr. Susama Patra
Additional Professor
Dept. of Pathology & Lab Medicine
AIIMS, Bhubaneswar-751019
HOD

Signature of Competent Authority.....

Date..... Designation of Officer.....